

## CM 20410 – CHLAMYDOSPORE AGAR

### INTENDED USE

For differentiation of *Candida albicans* from other *Candida* on the basis of chlamyospore formation.

### PRODUCT SUMMARY AND EXPLANATION

*Candida albicans* is a diploid sexual fungus (a form of yeast), and the causative agent of opportunistic oral and vaginal infections in humans. *C. albicans* is a commensal of skin, gastrointestinal and genitourinary tract. However, under certain conditions overgrowth of this results into oesopharyngeal candidiasis, vulvovaginal candidiasis and candidemia. Chlamyospores formation is the most differential characteristic of *C. albicans*. Chlamyospore Agar was specially designed for the differentiation of *C. albicans* from other species on the basis of chlamyospores formation. It is prepared according to the formula of Nickerson and Mankowsh.

**Test for chlamyospores:** Scratch cut mark like X onto the agar surface with inoculum using sterile needle. Aseptically place an alcohol-flamed and cooled cover slip onto the agar surface over the intersecting lines of the cut marks of X. Incubate plates at 20-25°C for 2-6 days. Temperature should not be higher than 25°C since it will not permit chlamyospore formation. Observe the plates under low power of microscope. After incubation, most strains of *C. albicans* and *C. stellatoide* will form typical chlamyospores. Chlamyospores will be seen along the edge of the cover slip. Chlamyospores are round, thick walled, blue coloured and at the terminal ends of hyphae. Some *C. albicans* strains may lose their ability to produce chlamyospores after repeated subculturing.

### COMPOSITION

Ingredients	Gms / Ltr
Ammonium sulphate	1.000
Monopotassium phosphate	1.000
Biotin	0.000005
Trypan blue	0.100
Purified polysaccharide	20.000
Agar	15.000

### PRINCIPLE

Ammonium sulphate acts as sources of ions that simulate metabolism. Monopotassium phosphate provides buffering to the medium. Biotin provides the necessary vitamins required for metabolism. Purified polysaccharide acts as a source of carbon. Trypan blue is a vital dye absorbed selectively by the chlamyospores and imparts blue colour to chlamyospores, whereas the filaments are colourless.

### INSTRUCTION FOR USE

- Dissolve 37.1grams in 1000 ml purified / distilled water.
- Heat to boiling to dissolve the medium completely.
- Sterilize by autoclaving at 15 psi pressure (121°C) for 15 minutes.
- Cool to 45-50°C. Mix well and pour into sterile Petri plates.

### QUALITY CONTROL SPECIFICATIONS



Appearance of Powder : Cream to blue homogeneous free flowing powder.  
 Appearance of prepared medium : Blue coloured opaque gel forms in Petri plates.  
 pH (at 25°C) : 5.1±0.2

#### INTERPRETATION

Cultural characteristics observed after incubation.

Microorganism	ATCC	Inoculum (CFU/ml)	Growth	Recovery	Chlamyospores	Incubation Temperature	Incubation Period
Candida albicans	10231	10-100	Good-luxuriant	>=50%	Positive	20-25°C	2-6 Days
Candida albicans	24408	10-100	Good-luxuriant	>=50%	Negative	20-25°C	2-6 Days
Candida tropicalis	1369	10-100	Good-luxuriant	>=50%	Negative	20-25°C	2-6 Days

#### PACKAGING:

In pack size of 100 gm and 500 gm bottles.

#### STORAGE

Dehydrated powder, hygroscopic in nature, store in a dry place, in tightly-sealed containers between 25-30°C and protect from direct sunlight. Under optimal conditions, the medium has a shelf life of 4 years. When the container is opened for the first time, note the time and date on the label space provided on the container. After the desired amount of medium has been taken out replace the cap tightly to protect from hydration.

Product Deterioration: Do not use if they show evidence of microbial contamination, discoloration, drying or any other signs of deterioration.













#### DISPOSAL

After use, prepared plates, specimen/sample containers and other contaminated materials must be sterilized before discarding.

#### REFERENCES

1. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.
2. Nickerson, 1953, J. Infect. Dis., 92:20.
3. Ryan K. J., Ray C. G., (Eds.), 2004, Sherris Medical Microbiology, 4th Ed., McGraw Hill.
4. Salfinger Y., and Tortorello M.L., 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.



 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative <small>MedNet GmbH Buckstraße 10 48163 Münster, Germany</small>	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.  
\*For LabUse Only

