

CM 23721 – BRUCELLA SELECTIVE SUPPLEMENT

INTENDED USE

For selective isolation of Brucella species from milk.

COMPOSITION

Ingredients	Concentration
Polymyxin B sulphate	2500IU
Bacitracin	12500IU
Nystatin	50000IU
Cycloheximide	50mg
Nalidixic acid	2.500mg
Vancomycin	10mg

(Per vial sufficient for 500 ml medium)

Applied in Medium

BLOODAGAR BASE No. 2

BRUCELLA AGAR BASE

BRUCELLA AGAR BASE (VEG.)

BRUCELLA AGAR BASE, MODIFIED

BRUCELLA BROTH BASE

BRUCELLA SELECTIVE MEDIUM BASE

BRUCELLA AGAR BASE W/ 1.0% DEXTROSE

COLUMBIA BLOOD AGAR BASE

COLUMBIA BLOOD AGAR BASE (VEG.)

BLOOD AGAR BASE NO. 2 (ISO 11290-1:2017)

BLOOD AGAR BASE NO. 2 WITH 1.2% AGAR

COLUMBIA BLOOD AGAR BASE W/ 1% AGAR













INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 10 ml of 50% methanol. Shake to form a uniform suspension. Add the contents to 500ml of sterile, molten, cooled (45-50°C) media such as, BLOOD AGAR BASE No. 2 / BLOOD AGAR BASE NO. 2 (ISO 11290-1:2017) / BLOOD AGAR BASE NO. 2 WITH 1.2% AGAR / COLUMBIA BLOOD AGAR BASE / COLUMBIA BLOOD AGAR BASE (VEG.) / COLUMBIA BLOOD AGAR BASE W/1% AGAR with 5-10% v/v inactivated horse serum and 1% w/v sterile dextrose / BRUCELLA AGAR BASE (VEG.) / BRUCELLA AGAR BASE, MODIFIED / BRUCELLA BROTH BASE / BRUCELLA SELECTIVE MEDIUM BASE / BRUCELLA AGAR BASE W/ 1.0% DEXTROSE with 5-10% v/v inactivated horse serum. Mix well and pour into sterile petri plates / tubes.

STORAGE

Vials should be stored in sealed container at 2-8°C.



 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative <small>Max-Mer GmbH Südstraße 10 81121 München, Germany</small>	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.
*For LabUse Only