

CM 23753 – CHROMOGENIC ENRICHMENT SUPPLEMENT

INTENDED USE

For the isolation, detection, and enumeration of *Listeria monocytogenes*.

COMPOSITION

Ingredients	Concentration
L – phosphatidylinositol	10ml

(Per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Aseptically add the contents of 2 vials of CHROMOGENIC ENRICHMENT SUPPLEMENT (TS 031) at a slightly elevated temperature (70°C) with vigorous shaking to 970ml of sterile molten CHROMOGENIC LISTERIA AGAR BASE / CHROMOGENIC L-MONO LISTERIA DIFFERENTIAL AGAR along with sterile rehydrated contents of 2 vials

of CHROMOGENIC LISTERIA SELECTIVE SUPPLEMENT (TS 205). Mix well and pour into sterile Petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MediStar GmbH Buckhorn 10, 48147 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.