

## **CM 22431 -SOYABEAN CASEIN DIGEST AGAR PLATE W/ 1% GLYCEROL,0.5% LECITHIN & 4% POLYSORBATE 80 (g - irradiated) (Triple Pack)**

### INTENDED USE

For determining of sanitization of containers,equipment, surfaces, water miscible cosmetics etc.

### PRODUCT SUMMARY AND EXPLANATION

Soyabean casein Digest Agar Plate w/ 1%Glycerol, 0.5%Lecithin and 4%Polysorbate 80 plates are recommended for the isolation of microorganisms from environmental surfaces and is used primarily to monitor microbial contamination, enumerate the number of microbial colonies growing on a variety of surfaces sanitized with quaternary ammonium compounds, and to assist in determining surface sanitation.

### COMPOSITION

Ingredients	Gms / Ltr
Polysorbate 80 (Tween 80)	40.000
Agar	15.000
Casein enzymic hydrolysate	15.000
Glycerol	10.000
Papaic digest of soyabean meal	5.000
Sodium chloride	5.000
Lecithin	5.000

### PRINCIPLE

Medium contains Casein enzymic hydrolysate and papaic digest of soyabean meal which helps to provide nitrogenous compounds and other nutrients essential for microbial replication. Sodium chloride is added to maintain cellular osmotic equilibrium. Lecithin and polysorbate 80 are added to the formulation to neutralize germicidal or disinfectant residues. Neutralization of these residues reduces their inhibitory effect which ultimately results in lowering of microbial count. Quaternary ammonia compounds are neutralized by lecithin, while phenolic disinfectants and hexachlorophene are neutralized by polysorbate 80. Together, lecithin and polysorbate 80 neutralize ethanol.

### INSTRUCTION FOR USE

Eitherstreak,inoculateorsurface spread the test inoculum aseptically on the plate. Alternatively, these plates can also be used as contact platesfor environmental monitoring.

### QUALITY CONTROL SPECIFICATIONS

Appearance	:	Light to medium amber colour medium.
Quantity of Medium	:	30±2 ml of medium in 90 mm sterile disposable Petri-plates.
pH (at 25°C)	:	7.3± 0.2
Dose of irradiation:	:	15-25 kGy
Sterility Check	:	Passes release criteria

### INTERPRETATION



Cultural characteristics observed after incubation.

Microorganism	ATCC	Inoculum (CFU/ml)	Growth	Recovery	Incubation Temperature	Incubation Period
Escherichia coli	25922	50-100	Luxuriant	>=70 %	35-37 °C	18-24 Hours
Staphylococcus aureus	25923	50-100	Luxuriant	>=70 %	35-37 °C	18-24 Hours
Pseudomonas aeruginosa	27853	50-100	Luxuriant	>=70 %	35-37 °C	18-24 Hours
Bacillus subtilis	6633	50-100	Luxuriant	>=70 %	35-37 °C	18-24 Hours
Salmonella typhimurium	14028	50-100	Luxuriant	>=70 %	35-37 °C	18-24 Hours
Candida albicans	10231	50-100	Luxuriant	>=70 %	30-35 °C	48-72 Hours

#### PACKAGING:

Triplelayeredpacking containing 5 No. of plates with one silica gel desiccant bag packed inside it.

#### STORAGE

Onreceipt, store the plates at 15–30 °C. Avoid freezing and overheating. Do not open until ready to use. Prepared plates stored in their original sleeve wrapping until just prior to use may be inoculated up to the expiration date and incubated for recommended incubation times. Allow the medium to warm to room temperature before inoculation. Product Deterioration: Do not use plates if they show evidence of microbial contamination, discoloration, drying, cracking or other signs of deterioration.

#### DISPOSAL

Afteruse,prepared plates, specimen/sample containers and other contaminated materials must be sterilized before discarding.

#### REFERENCES

- TheUnited States Pharmacopoeia. 2009. Amended Chapters 61, 62 & 111, The United States Pharmacopoeial Convention Inc., Rockville, MD.
- Directorate for the Quality of Medicines of the Council of Europe (EDQM). 2009. The European Pharmacopoeia, Amended Chapters 2.6.12, 2.6.13, 5.1.4, Council of Europe, 67075 Strasbourg Cedex, France.
- Japanese Pharmacopoeia. 2008. Society of Japanese Pharmacopoeia. Amended Chapters 35.1, 35.2.7. The Minister of Health, Labor, and Welfare.
- Indian Pharmacopoeia. 2010. Govt. of India, Ministry of Health and Family Welfare, New Delhi, India.
- MacFaddin, J.F. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. I. Williams & Wilkins, Baltimore.
- Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis.
- Chapin, K.C., and P.R. Murray. 1999. Media, p. 1687-1707. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 7th ed. American Society for Microbiology, Washington, D.C
- Clesceri, L.S., A.E. Greenberg, and A.D. Eaton (ed.). 1998. Standard methods for the examination of water and wastewater, 20th ed. American Public Health Association, Washington, D.C.
- Downes, F.P. and K. Ito. (ed.). 2001. Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington, D.C.
- ISO 11137-1: 2006 + Amd 1:2013. Sterilization of health care products – Radiation - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices.
- ISO 11137-2:2013. Sterilization of health care products -- Radiation -- Part 2: Establishing the sterilization dose.



Quantity



Lot / Batch Number



Temperature Unit



Manufacturer



Best Before



Certification of Good Manufacturing Practices



Catalogue No.



Authorized Representative



European Conformity



QR Code



Consults Instructions for use :



For In Vitro Diagnostic Use

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.



\*For Lab Use Only

th

