

CM 23712 – BACILLUS CEREUS SELECTIVE SUPPLEMENT

INTENDED USE

For selective isolation of Bacillus species.

COMPOSITION

Ingredients	Concentration
Polymixin B	50000IU

(Per vial sufficient for 500 ml medium)

Applied in Medium


CHROMOGENIC BACILLUS CEREUS AGAR

INSTRUCTION FOR USE

Aseptically reconstitute 1 vial with 5 ml of sterile distilled water. Mix well and aseptically add it to 500 ml sterile, molten CHROMOGENIC BACILLUS CEREUS AGAR. Mix well and pour into sterile petri plates / tubes.

STORAGE

Vials should be stored in sealed container at 2-8°C.



NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For LabUse Only

