

CM 20791 - HARTLEY'S DIGEST BROTH (VEG.)

INTENDED USE

For cultivation of various bacteria from blood especially Streptococci and Corynebacterium diphtheriae.

PRODUCT SUMMARY AND EXPLANATION

Hartley's Digest Broth (Veg) is prepared by using Veg special infusion in place of tryptic digest of heart muscle thus making the medium free of BSE/TSE risks. This is the modification of the medium described by Hartley for the production of diphtheria toxin. Hartley's Digest Broth (Veg) can be used for the recovery of fastidious organisms such as Group A, C, G Streptococci and Streptococcus pneumoniae from small inocula. It can also be used to recover Corynebacterium diphtheriae. Monckton used it as enrichment medium for the same organism. It may be used for cultivation of blood samples, sterility testing, production of diphtheria toxin etc.

COMPOSITION

Ingredients	Gms / Ltr
Veg special infusion	29.000

PRINCIPLE

Hartley's Digest Broth (Veg) is prepared as per the method described by Cruickshank. Veg special infusion provides nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other essential nutrients.

INSTRUCTION FOR USE

Dissolve 29.0 grams in 1000 ml distilled water.

Heat if necessary to dissolve the medium completely.

Dispense in tubes or flasks as desired and sterilize by autoclaving at 15 psi pressure (121°C) for 15 minutes.

QUALITY CONTROL SPECIFICATIONS

Appearance of Powder	: Light yellow coloured, may have slightly greenish tinge, homogeneous, free flowing powder.
Appearance of prepared medium	: Light yellow coloured, clear solution without any precipitate.
pH (at 25°C)	: 7.6±0.2

INTERPRETATION

Cultural characteristics observed after an incubation.

Microorganism	ATCC	Inoculum (CFU/ml)	Growth	Incubation Temperature	Incubation Period
Corynebacterium diphtheriae	11913	50-100	Luxuriant	35 - 37°C	24-48 Hours
Enterococcus faecalis	29212	50-100	Luxuriant	35 - 37°C	24-48 Hours



Staphylococcus aureus	25923	50-100	Luxuriant	35 - 37°C	24-48 Hours
Streptococcus pneumoniae	6303	50-100	Luxuriant	35 - 37°C	24-48 Hours
Streptococcus pyogenes	19615	50-100	Luxuriant	35 - 37°C	24-48 Hours

PACKAGING:

In pack size of 100 gm and 500 gm bottles.

STORAGE

Dehydrated powder, hygroscopic in nature, store in a dry place, in tightly-sealed containers between 25-30°C and protect from direct sunlight. Under optimal conditions, the medium has a shelf life of 4 years. When the container is opened for the first time, note the time and date on the label space provided on the container. After the desired amount of medium has been taken out replace the cap tightly to protect from hydration.

Product Deterioration: Do not use if they show evidence of microbial contamination, discoloration, drying or any other signs of deterioration.

DISPOSAL

After use, prepared plates, specimen/sample containers and other contaminated materials must be sterilized before discarding.

REFERENCES

1. Hartley P., 1922, J. Path. Bact., 25:479.
2. Monckton J.C., 1947, Bull. Inst. Med. Lab. Technol., 13(1):2.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MediMet GmbH Buckstrasse 10 48143 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only

