

CM 23725 – CAMPYLOBACTER SELECTIVE SUPPLEMENT (Karmali)

INTENDED USE

For isolating thermotolerant Campylobacter species.

COMPOSITION

Ingredients	Concentration
Sodium pyruvate	50mg
Vancomycin	10mg
Cefoperazone	16mg
Cycloheximide	50mg

(Per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 5.0 ml of 50% ethanol. Mix well and aseptically add to 500 ml sterile, molten, cooled (45-50°C) CHROMOGENIC CAMPYLOBACTER AGAR BASE or to KARMALI CAMPYLOBACTER AGAR BASE along with 5 ml of Hemin Solution (16 mg / 5ml). Mix well and pour in sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MediNet GmbH Buckstrasse 10, 48143 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only

