

CM 23727 – CAMPYLOBACTER SELECTIVE SUPPLEMENT (With Hemin) (Karmali)

INTENDED USE

For isolating thermotolerant Campylobacter species.

COMPOSITION

Ingredients	Concentration
Sodium pyruvate	50mg
Vancomycin	10mg
Cefoperazone	16mg
Cycloheximide	50mg
Hemin	16mg

(Per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 5.0 ml of 95% ethanol. Mix well and aseptically add to 500 ml sterile, molten, cooled (45 - 50°C) KARMALI CAMPYLOBACTER AGAR BASE. Mix well and pour in sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative Müller Martini GmbH Bockhorn 10 61163 Frankfurt, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only

