

CM 23729 – CAMPYLOBACTER SUPPLEMENT-IV (Preston)

INTENDED USE

For isolation of Campylobacter species.

COMPOSITION

Ingredients	Concentration
Polymyxin B sulphate	2500IU
Rifampicin	5mg
Trimethoprim	5mg
Cycloheximide	50mg

(Per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% aqueous acetone. Mix well and aseptically add it to 475 ml of sterile, molten, cooled (45-50°C) TM 603 - CAMPYLOBACTER ENRICHMENT BROTH BASE (PRESTON ENRICHMENT BROTH BASE) / CAMPYLOBACTER ENRICHMENT AGAR BASE (PRESTON ENRICHMENT AGAR BASE) / PRESTON AGAR BASE / PRESTON AGAR BASE (VEG.) along with 25 ml sterile, lysed horse blood. Mix well and dispense into sterile tubes / petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative Mallatier GmbH Buckhorn 10, 40147 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only