

CM 23732 – CAMPYLOBACTER SUPPLEMENT-II (Butzler)

INTENDED USE

For selective isolation of Campylobacter species.

COMPOSITION

Ingredients	Concentration
Bacitracin	12500IU
Colistin sulphate	5000Units
Cycloheximide	25mg
Cephazolin sodium	7.500mg
Novobiocin	2.500mg

(Per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 5 ml of 50% aqueous ethanol and mix well to dissolve. Avoid frothing of the solution. Aseptically add the contents to 500 ml of sterile, molten, cooled (45-50°C) COLUMBIA BLOOD AGAR BASE / COLUMBIA BLOOD AGAR BASE (VEG.) / COLUMBIA BLOOD AGAR BASE W/ 1% AGAR / BLOOD AGAR BASE No. 2 / BLOOD AGAR BASE NO. 2 (ISO 11290-1:2017) / BLOOD AGAR BASE NO. 2 WITH 1.2% AGAR / BRUCELLA AGAR BASE / BRUCELLA AGAR BASE (VEG.) along with 5% v/v defibrinated Sheep blood or 5-7% v/v laked horse blood and rehydrated contents of 1 vial of Campylobacter Growth Supplement. Mix gently and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.