

CM 23733 – CAMPYLOBACTER SUPPLEMENT-III (Skirrow)

INTENDED USE

For selective isolation of Campylobacter species.

COMPOSITION

Ingredients	Concentration
Vancomycin	5mg
Polymyxin B sulphate	1250IU
Trimethoprim	2.500mg

(Per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 2 ml sterile distilled water and mix well to dissolve. Avoid frothing of the solution. Aseptically add the rehydrated contents to 500 ml of sterile, molten, cooled (45-50°C) COLUMBIA BLOOD AGAR BASE / COLUMBIA BLOOD AGAR BASE (VEG.) / COLUMBIA BLOOD AGAR BASE W/ 1% AGAR / BLOOD AGAR BASE No. 2 / BLOOD AGAR BASE NO. 2 (ISO 11290-1:2017) / BLOOD AGAR BASE NO. 2 WITH 1.2% AGAR / CAMPYLOBACTER AGAR BASE with 5-7% v/v lysed defibrinated horse blood BRUCELLA AGAR BASE / BRUCELLA AGAR BASE (VEG.) / BRUCELLA AGAR BASE, MODIFIED with 5% inactivated horse serum. Mix gently and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative Müller Gekit Sackstrasse 10, 41543 Alsdorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.