

CM 23736 – CATC SUPPLEMENT

INTENDED USE

Recommended for the selective isolation and detection of *Enterococcus faecalis* by means of TTC reduction.

COMPOSITION

Ingredients	Concentration
Sodium carbonate	1g
2,3,5 Triphenyltetrazolium chloride	0.050g
Sodium azide	0.200g

(Per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 10 ml sterile distilled water and add it to 500 ml sterile, molten, cooled (45-50°C) CITRATE AZIDE TWEEN CARBONATE BASE. Mix well and pour into sterile Petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative Muller GmbH Buckhorn 10, 40112 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only