

CM 23737 – CC SUPPLEMENT

INTENDED USE

Recommended for the selective cultivation of pathogenic fungi.

COMPOSITION

Ingredients	Concentration
Cycloheximide	200mg
Chloramphenicol	25mg

(Per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 5 ml of distilled water. Mix well and aseptically add it to 500 ml sterile, molten, cooled (45-50°C) DEXTROSE AGAR BASE, EMMONS (SABOURAUD DEXTROSE AGAR BASE, MODIFIED). Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedWet GmbH Seidenstrasse 15, 40147 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only

