

CM 23739 – CEFIXIME SUPPLEMENT

INTENDED USE

Recommended for selective isolation of E.coli O157:H7.

COMPOSITION

Ingredients	Concentration
Cefixime	0.025mg

(Per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial with 2 ml sterile distilled water. Mix well and aseptically add one vial to 500 ml or two vials to 1000 ml of sterile, molten cooled (45-50°C) - MacCONKEY SORBITOL AGAR BASE W/ RHAMNOSE as desired. Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWet GmbH Seckestraße 116 48163 Münster, Germany	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only