

CM 23743 – CHLORAMPHENICAL SELECTIVE SUPPLEMENT

INTENDED USE

For the selective isolation of yeasts and moulds.

COMPOSITION

Ingredients	Concentration
Chloramphenicol	50mg

(Per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the content of 1 vial aseptically with 2 ml of ethanol. Mix well and aseptically add to 500 ml of sterile molten, cooled (45-50°C) SABHI AGAR BASE / ASPERGILLUS DIFFERENTIATION MEDIUM BASE / DICHLORAN GLYCEROL MEDIUM BASE W/ROSE BENGAL / DICHLORAN GLYCEROL MEDIUM BASE W/ROSE BENGAL (VEG.) / ROSE BENGAL AGAR BASE. Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 Best Before	 QTY. Quantity	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 LOT/ B. NO. Lot / Batch Number	 Consults Instructions for Use	 QR Code	

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only

