

CM 23747 – CHROMOGENIC EC O157:H7 SELECTIVE SUPPLEMENT, MODIFIED

INTENDED USE

For selective isolation and detection of *Enterococcus faecium*.

COMPOSITION

Ingredients	Concentration
Cephalexin	25mg
Aztreonam	37.500mg

(Per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and add aseptically to 500 ml sterile, molten, cooled (45-50°C) CHROMOGENIC ENTEROCOCCUS FAECIUM AGAR BASE. Mix well and pour into sterile Petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MediStar GmbH Buckhorn 16, 48143 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only