

CM 23748 – CHROMOGENIC CANDIDA DIFFERENTIAL SELECTIVE SUPPLEMENT

INTENDED USE

Recommended for rapid and direct isolation and identification of Candida species from mixed cultures.

COMPOSITION

Ingredients	Concentration
Chloramphenicol	250mg

(Per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 2 ml of 95% ethanol. Mix well and aseptically add to 500 ml of sterile, molten cooled (45-50°C) CHROMOGENIC CANDIDA DIFFERENTIAL AGAR BASE. Mix well and pour into sterile Petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedStar GmbH Buckhorn 10, 40162 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only