

CM 23752 – CHROMOGENIC ECC SELECTIVE SUPPLEMENT

INTENDED USE

For the selective isolation of Escherichia coli and coliforms from water and food samples.

COMPOSITION

Ingredients	Concentration
Cefsulodin	10mg

(Per vial sufficient for 1000ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 2 ml sterile distilled water. Mix well and aseptically add to 1000 ml of sterile, molten, cooled (45-50°C) CHROMOGENIC ECC SELECTIVE AGAR BASE. Mix well and pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWet GmbH Seckestraße 11 48153 Münster, Germany	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only