

CM 23763 – COLISTIN SELECTIVE SUPPLEMENT

INTENDED USE

For selective isolation of Clostridium perfringens.

COMPOSITION

Ingredients	Concentration
Colistin sulphate	13 mg/L

(per vial sufficient for 1000 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water. Mix well and aseptically add it to 1000 ml of sterile, molten, cooled (45-50°C) CPC AGAR BASE W/ 1% CELLOBIOSE (TM 2027). Mix well and pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedNet GmbH Buckstraße 10 48163 Münster, Germany	 CE European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only