

CM 23772 – EGG YOLK EMULSION (100 ml/vl)

INTENDED USE

For use in various culture media.

COMPOSITION

Ingredients	Concentration
Egg yolk	30ml
Sterile saline	70ml

(per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Warmup (at 40-45°C) the refrigerated Egg Yolk Tellurite Emulsion with gentle shaking, to get a uniform emulsion. Aseptically add 50 ml of this suspension with 950 ml of sterilized, cooled (at 45-50°C), molten media BACILLUS CEREUS AGAR BASE / BAIRD PARKER AGAR BASE / BAIRD PARKER AGAR BASE(VEG.) / BAIRD PARKER AGAR BASE (RPF) (ISO 6888-1 & 2:1999) / BACILLUS CEREUS SELECTIVE AGAR BASE (MYP) (ISO 7932:2004) / BAIRD PARKER AGAR BASE (IS : 5887 (Part II) 1976, reaffirmed 2005.) / BAIRD PARKER AGAR BASE W/ SULPHA / C. BOTULINUM ISOLATION AGAR BASE / CLOSTRIDIUM PERFRINGENS AGAR / EGG YOLK AGAR BASE / EGG YOLK AGAR BASE, MODIFIED / MANNITOL EGG YOLK POLYMYXIN AGAR (ISO 7932:2004, ISO 21871:2006) / TRYPTOSE SULPHITE CYCLOSERINE AGAR BASE / PERFRINGENS AGAR BASE (T.S.C./S.F.P. AGAR BASE) (VEG.) / MYP AGAR BASE (PHENOL RED EGG YOLK POLYMYXIN AGAR BASE) (IS : 5887 (Part VI) 1976, reaffirmed 2005) / MYP AGAR BASE (PHENOL RED EGG YOLK POLYMYXIN AGAR BASE) (VEG.) / MANNITOL SALT AGAR BASE / MANNITOL SALT AGAR BASE (VEG.) / K.R.A.N.E.P. AGAR BASE / L.D. EGG YOLK AGAR BASE / KG AGAR BASE / LIPOVITELLIN SALT MANNITOL AGAR BASE / MCCLUNG TOABE AGAR BASE / TPEY AGAR BASE / ANAEROBIC EGG AGAR BASE / POLYMYXIN PYRUVATE EGG YOLK MANNITOL BROMOTHYMOL BLUE AGAR BASE (PEMBA)) / WILLIS AND HOBBS MEDIUM BASE (IS : 5887 (Part 4) : 1999) / MANNITOL SALT AGAR BASE (as per USP/BP/EP/JP) / MANNITOL SALT AGAR BASE (as per USP/BP/EP/JP) (VEG.). Mix well and dispense as desired.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedWor GmbH Seidenstr. 10 40152 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

