

CM 23774 – EGG YOLK TELLURITE EMULSION (100 ml/vl)

INTENDED USE

For identification of Staphylococci.

COMPOSITION

Ingredients	Concentration
Egg yolk	30ml
Sterile saline	64ml
Sterile 3.5% potassium tellurite solution	6ml

(per vial sufficient for 100 ml medium)

INSTRUCTION FOR USE

Warm up the refrigerated Egg Yolk Tellurite Emulsion to 40-45°C. Shake well to attain uniform emulsion (since on refrigeration, the emulsion has a tendency to form layers or small lumps). Aseptically add 50 ml in 950 ml of sterile, molten, cooled (45-50°C) Baird parker agar base (RPF) (ISO 6888-1 & 2:1999) / Baird parker agar base (IS : 5887 (Part II) 1976, reaffirmed 2005.) / Baird parker agar base (agar medium o) (as per USP/EP) / Baird parker agar base/ sulpha / Baird parker agar base / Baird parker agar base (veg.) / Chromogenic staphylococcus aureus agar base / Baird parker agar base (as per EP/IP/BP). Mix well and pour

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

