

CM 23779 – FRASER SELECTIVE SUPPLEMENT

INTENDED USE

For selective isolation and enrichment of *Listeria monocytogenes* from food and environmental samples.

COMPOSITION

Ingredients	Concentration
Acriflavine hydrochloride	25mg
Nalidixic acid	20mg

(per vial sufficient for 1000 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water. Mix well and aseptically add to 990 ml of sterile, cooled (45-50°C) FRASER SECONDARY ENRICHMENT BROTH BASE / FRASER BROTH BASE (ISO 11290-1: 2017) HALF FRASER BROTH BASE (ISO 11290-1:2017) / FRASER BROTH BASE. Mix well and dispense as desired.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedMer GmbH Südstraße 15 48143 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.