

CM 23781 – FRIIS SUPPLEMENT

INTENDED USE

A selective supplement for the detection of non-avian Mycoplasmas in pharmaceutical products in accordance with the European Pharmacopoeia.

COMPOSITION

Ingredients	Concentration
Bacitracin	125.0mg/vl
Methicillin	125.0mg/vl
Nalidixic acid	10.0mg/vl

(per vial sufficient for 1000ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 2 ml of 0.2N NaOH and 8 ml of sterile distilled water. Mix well and aseptically add to 1000 ml of sterile, cooled (45-50°C) FRIIS LIQUID MEDIUM BASE. Mix well and dispense as desired.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedWet GmbH Bachstrasse 15, 40142 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only

