

CM 23794 – HEAMOGLOBIN POWDER

INTENDED USE

For isolation of Neisseria spp.

COMPOSITION

Ingredients	Concentration
Haemoglobin powder	50G/ 100G

INSTRUCTION FOR USE

Aspecially prepared powder whose 2% w/v aqueous solution is autoclavable. The aqueous solution is chocolate brown, opaque and contains flocculent dispersible precipitate. It is used for 500 ml medium preparation of GC AGAR BASE - 5 gms / GC AGAR BASE (VEG.) - 5 gms / THAYER MARTIN MEDIUM BASE – 5 gms / CHOCOLATE NO. 2 AGAR BASE – 5 gms / TELLURITE BLOOD AGAR BASE – 10 gms / CHOCOLATE AGAR BASE - 10 gms / TRANSGROW MEDIUM BASE - 2 gms.

STORAGE

Vialsshouldbe stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative Muller GmbH Ruckstrasse 16 40143 Essen, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For LabUse Only

