

CM 23798 – KANAMYCIN SULPHATE SELECTIVE SUPPLEMENT

INTENDED USE

For selective isolation of Enterococci.

COMPOSITION

Ingredients	Concentration
Kanamycin sulphate	10mg

(per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C) KANAMYCIN ESCULIN AZIDE AGAR BASE / KANAMYCIN ESCULIN AZIDE BROTH BASE. Mix well and aseptically pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWet GmbH Sanktstraße 115 48143 Münster, Germany	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only