

CM 23801 – L. MONO ENRICHMENT SUPPLEMENT I

INTENDED USE

For selective differentiation of *Listeria monocytogenes* from other *Listeria* species.

COMPOSITION

Ingredients	Concentration
L-phosphatidylinositol	1.0 g
Distilled water	25 ml

(per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Thaw the contents of 1 vial of L. mono Enrichment Supplement I at room temperature. Aseptically add sterile contents of 1 vial of L. mono Enrichment Supplement I to 460 ml of sterile, molten L. mono Differential Agar Base along with sterile rehydrated contents of L. Mono Selective Supplement I (TS 227) and L. mono Selective Supplement II Mixwell and pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedNet GmbH Barkstrasse 10, 48163 Münster, Germany	 CE European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only