

CM 23808 – LANDER AF SUPPLEMENT B

INTENDED USE

For selective enrichment and transport of *Campylobacter* spp.

COMPOSITION

Ingredients	Concentration
Cyclohexamide	100mg/L

(Per vial sufficient for 1000ml medium)

INSTRUCTION FOR USE

Rehydrate the contents aseptically of 1 vial into to 1000 ml of sterile, molten LANDER TRANSPORT & STORAGE MEDIUM along with Lander VPT supplement A. Mix carefully and dispense as desired.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWet GmbH Seestraße 11 40113 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only