

CM 23811 – LECITHIN SOLUTION

INTENDED USE

For based on PCPLC activity.

COMPOSITION

Ingredients	Concentration
Soya lecithin	20ml

(Per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Thaw the contents of one vial at room temperature. Mix well and aseptically add it to 480 ml of sterile molten, cooled (45-50°C) CHROMOGENIC L. MONO DIFFERENTIAL AGAR BASE along with sterile contents of 1 vial of Modified L. Mono Selective Supplement. Mix well and pour into sterile Petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWet GmbH Seidenstrasse 15, 40142 Düsseldorf, Germany	 CE European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only