

CM 23812 – LEGIONELLA (GVPC) SELECTIVE SUPPLEMENT

INTENDED USE

For isolation of Legionella species.

COMPOSITION

Ingredients	Concentration
Glycine	1.500g
Vancomycin hydrochloride	0.500mg
Polymyxin B sulphate	40000IU
Cycloheximide	40mg

(Per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water. Mix well and aseptically add to 440 ml sterile, molten, cooled (45-50°C) LEGIONELLA AGAR BASE / LEGIONELLA AGAR BASE W/O CHARCOAL along with 50 ml of Legionella Growth Supplement BCYE or to 500 ml sterile, molten, cooled (45-50°C) BUFFERED CHARCOAL YEAST EXTRACT AGAR MEDIUM (ISO 11731-2:2017). Mix well and pour into sterile Petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWet GmbH Bachstrasse 16, 48147 Münster, Germany	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For LabUse Only

