

CM 23820 – LEGIONELLA SELECTIVE SUPPLEMENT-III

INTENDED USE

For selective isolation of Legionella species from mixed cultures.

COMPOSITION

Ingredients	Concentration
Polymyxin B sulphate	50000Unit
Glycine	1.500mg
Vancomycin	2.500mg

(Per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 500 ml of sterile, molten, cooled (45-50°C) BUFFERED CHARCOAL YEAST EXTRACT AGAR BASE along with rehydrated contents of Legionella Selective Supplement IV and Legionella Supplement. Mix well and pour into sterile petri plates.

STORAGE

Vialsshouldbe stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative Mallatier GmbH Buckhorn 16, 40142 Mönchen, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For LabUse Only