

CM 23826 – LISTERIA MOXALACTAM SUPPLEMENT MODIFIED

INTENDED USE

Recommended for the selective isolation of Listeria species from food samples in accordance with FDA BAM,1998.

COMPOSITION

Ingredients	Concentration
Colistin sulfate	5mg
Moxalactum	10mg

(Per vial sufficient for 500 ml medium)

Applied in Medium

TM 2162 - LISTERIA OXFORD AGAR BASE W/ 1.2% AGAR

TM 2163 - LISTERIA OXFORD MEDIUM BASE, MODIFIED

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial with 5 ml of sterile distilled water and aseptically to 500 ml of sterile molten cooled (45-50°C) LISTERIA OXFORD AGAR BASE W/ 1.2% AGAR / LISTERIA OXFORD MEDIUM BASE, MODIFIED. Mix well and pour into sterile Petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative Mueller GmbH Sackstrasse 10 48143 Münster, Germany	 European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For LabUse Only