

CM 23832 – M-CP SELECTIVE SUPPLEMENT I

INTENDED USE

Recommended for selective isolation of Clostridium perfringens.

COMPOSITION

Ingredients	Concentration
D-Cycloserine	200mg
Polymyxin B sulphate	12.500mg

(Per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to 485 ml sterile, molten, cooled (45-50°C) M-CP AGAR BASE along with one vial of M-CP Selective Supplement II / M-CP Selective Supplement, Modified. Mix well and pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedWet GmbH Buckhorn 16, 48163 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only