

CM 23833 – M-CP SELECTIVE SUPPLEMENT II

INTENDED USE

Recommended for selective isolation of Clostridium perfringens.

COMPOSITION

Ingredients	Concentration
Phenolphthalein diphosphate	0.050g
Distilled water	10ml

(Per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Warm up the refrigerated 0.5% phenolphthalein diphosphate solution to room temperature and add aseptically 10 ml of solution to 485 ml sterile, molten, cooled (45-50°C) M-CP AGAR BASE along with rehydrated contents of one vial of M-CP Selective Supplement I / M-CP Selective Supplement, Modified. Mix well and pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MediStar GmbH Buckenhof 16, 48147 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only