

## CM 23834 – M-CP SELECTIVE SUPPLEMENT, MODIFIED

### INTENDED USE

Recommended for selective isolation of Clostridium perfringens.

### COMPOSITION

| Ingredients                 | Concentration |
|-----------------------------|---------------|
| Phenolphthalein diphosphate | 0.050g        |

(Per vial sufficient for 500ml medium)

### INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water. Mix well. Add aseptically to 485 ml sterile, molten, cooled to (45-50°C) M-CP AGAR BASE along with rehydrated contents of one vial of M-CP Selective Supplement I and one vial of M-CP Selective Supplement II. Mix well and pour into sterile petri plates.

### STORAGE

Vial should be stored in sealed container at 2-8°C.

|                                                                                                                                         |                                                                                                                                                                                                             |                                                                                                            |                                                                                                                             |                                                                                                                           |                                                                                                       |
|-----------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------|
| <br>GMP<br>Good Manufacturing<br>Practices Certified | <br>IVD<br>For In Vitro Diagnostic Use                                                                                   | <br>QTY.<br>Quantity    | <br>LOT/<br>B. NO.<br>Lot / Batch Number | <br>REF<br>Catalogue Number          | <br>Manufacturer |
| <br>Temperature Unit                                 | <br>EC REP<br>Authorized Representative<br><small>MedWet GmbH<br/>Steinbockstr. 11<br/>40113 Düsseldorf, Germany</small> | <br>European Conformity | <br>QR<br>Code                           | <br>Consults instructions<br>for Use | <br>Best Before  |

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

\*For Lab Use Only