

CM 23835 – M-PA SELECTIVE SUPPLEMENT

INTENDED USE

Recommended for selective isolation of *Pseudomonas aeruginosa*.

COMPOSITION

Ingredients	Concentration
Sulphapyridine	176mg
Kanamycin	8.500mg
Nalidixic acid	37mg
Amphotericin B	5mg

(Per vial sufficient for 1000ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to 1000 ml of sterile, molten-cooled (45-50°C) M-PA AGAR BASE. Mix well and aseptically pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative Mallatier GmbH Buckhorn 10, 40143 Essen, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only