

CM 23836 – MCBRIDE LISTERIA SUPPLEMENT

INTENDED USE

For selective isolation of *Listeria* species.

COMPOSITION

Ingredients	Concentration
Cycloheximide	200mg

(Per vial sufficient for 1000 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial with 10 ml sterile distilled water. Mix well and aseptically add to 990 ml of sterile, molten, cooled (45-50°C) MCBRIDE LISTERIA AGAR BASE along with 5% v/v sterile defibrinated blood. Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative Muehler GmbH Seckestraße 116 48153 Münster, Germany	 CE European Conformity	 QR Code	 Consults instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only