

CM 23842 – MODIFIED CFC SELECTIVE SUPPLEMENT

INTENDED USE

For selective isolation of Pseudomonas species.

COMPOSITION

Ingredients	Concentration
Cetrimide	5mg
Sodium fusidate	5mg
Cephalothin sodium salt	25mg

(Per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 5 ml distilled water. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C), CFC AGAR BASE (CEPHALOTHIN-SODIUM FUSIDATE-CETRIMIDE AGAR) / CFC AGAR BASE (CEPHALOTHIN-SODIUM FUSIDATE-CETRIMIDE AGAR). Mix well and pour into sterile Petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWet GmbH Sankt Leon 11 40113 Alsdorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only