

CM 23844 – MODIFIED L.MONO SELECTIVE SUPPLEMENT

INTENDED USE

For selective isolation of L.monocytogenes and Listeria species based on PCPLC activity.

COMPOSITION

Ingredients	Concentration
Nalidixic acid	0.10gm
Polymyxin B Sulfate	0.000004564gm
Ceftazidime	0.10gm
Amphotericin	0.005gm

(Per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to 480 ml of sterile molten, cooled (45-50°C) CHROMOGENIC L. MONO DIFFERENTIAL AGAR BASE along with sterile contents of 1 vial of Lecithin solution. Mix well and pour into sterile Petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative Mallatier GmbH Buckhorn 10, 48143 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only