

CM 23858 – NOVOBIOCIN SUPPLEMENT

INTENDED USE

For presumptive identification of *Aeromonas hydrophila*.

COMPOSITION

Ingredients	Concentration
Novobiocin	5mg

(per vial sufficient for 1000ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water and aseptically add to 1000 ml of sterile, molten, cooled (45-50°C) MULLER KAUFFMAN TETRATHIONATE BROTH BASE / SEMISOLID RAPPAPORT VASSILIADIS MEDIUM, MODIFIED / SOYABEAN BILE BROTH BASE / DOUBLE MODIFIED LYSINE IRON AGAR BASE / RS MEDIUM BASE. Mix well and pour into sterile Petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative Muller GmbH Buckhorn 16 48143 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only