

CM 23860 – OLEIC ALBUMIN SUPPLEMENT

INTENDED USE

For cultivation of Mycobacterium tuberculosis.

COMPOSITION

Ingredients	Concentration
Alkalinized oleic acid	0.050%
Albumin fraction V	5%

(per vial sufficient for 200 ml medium)

INSTRUCTION FOR USE

Warm up the refrigerated supplement to room temperature. Shake well and aseptically add the entire content (20 ml) of 1 vial of sterile Oleic Albumin Supplement to 180 ml of sterile, molten, cooled (45-50°C) DUBOS OLEIC BROTH BASE / DUBOS OLEIC AGAR BASE along with 5000 to 10000 units of penicillin. Mix well.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative Müller-Merz GmbH Buckhorn 10 48143 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only