

CM 23862 – OXFORD LISTERIA SUPPLEMENT

INTENDED USE

For selective isolation of Listeria species.

COMPOSITION

Ingredients	Concentration
Cycloheximide	200mg
Colistin sulphate	10mg
Acriflavin	2.500mg
Cefotetan	1mg
Fosfomycin	5mg

(per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix gently to dissolve and aseptically add to 500 ml of sterile, molten, cooled (45-50°C) LISTERIA OXFORD MEDIUM BASE (ISO 11290-1) / LISTERIA OXFORD MEDIUM BASE, MODIFIED. Mix well and pour into sterile Petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative Muller GmbH Buckhorn 10, 40143 Essen, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only