

CM 23865 – PARK AND SANDERS SELECTIVE SUPPLEMENT B

INTENDED USE

For selective isolation and enumeration of thermotolerant Campylobacter species.

COMPOSITION

Ingredients	Concentration
Cefoperazone	32mg
Cycloheximide	100mg

(per vial sufficient for 950ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 5 ml of 50% aqueous ethanol. Mix well and aseptically add it to 950 ml of sterile, cooled (45-50°C) CAMPYLOBACTER CEFEX BROTH BASE / PARK & SANDER ENRICHMENT BROTH BASE / PARK & SANDER ENRICHMENT BROTH BASE (ISO/DIS 10272:1995) along with 50 ml of sterile defibrinated lysed horse blood and rehydrated contents of 1 vial of Park and Sanders Selective Supplement A / CAMPYLOBACTER CEFEX AGAR BASE containing 50 ml sterile defibrinated sheep blood. Mix well and dispense.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWet GmbH Seckestraße 15 40162 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only