

CM 23868 – S.F. P. SUPPLEMENT (Perfringens S.F.P. Supplement)

INTENDED USE

For selective isolation of Clostridium perfringens.

COMPOSITION

Ingredients	Concentration
Kanamycin sulphate	6mg
Polymyxin B	15000IU

(per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 2 ml sterile distilled water. Mix well and aseptically add it to 475 ml of sterile, molten, cooled (45-50°C) PERFRINGENS AGAR BASE (T.S.C./S.F.P. AGAR BASE) / PERFRINGENS AGAR BASE (T.S.C./S.F.P. AGAR BASE) (VEG.) / TRYPTOSE SULPHITE CYCLOSERINE AGAR BASE / S.F.P. AGAR BASE along with 25 ml of sterile Egg Yolk Emulsion. Mix well and pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative MedWor GmbH Buckhorn 16 40162 Düsseldorf, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only