

CM 23869 – PERFRINGEN'S T.S.C. SUPPLEMENT (TSC SUPPLEMENT)

INTENDED USE

For the selective isolation of Clostridium perfringens.

COMPOSITION

Ingredients	Concentration
D-Cycloserine	200mg

(per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 2 ml sterile distilled water. Mix well and aseptically add to 475 ml of sterile, molten, cooled (45-50°C) PERFRINGENS AGAR BASE (T.S.C./S.F.P. AGAR BASE) / PERFRINGENS AGAR BASE (T.S.C./S.F.P. AGAR BASE) (VEG.) / S.F.P. AGAR BASE along with 25 ml of Egg Yolk Emulsion or 500 ml of sterile, molten, cooled (45-50°C) PERFRINGENS AGAR BASE (TRYPTOSE SULPHITE CYCLOSERINE AGAR BASE) (ISO 7937: 2004, ISO 14189:2013) / TRYPTOSE CYCLOSERINE AZIDE AGAR BASE / TRYPTOSE CYCLOSERINE DEXTROSE AGAR BASE / TRYPTONE YEAST SODIUM SULPHITE AGAR BASE (ISO 14189:2013). Mix well and pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedWet GmbH Seckestraße 15 48161 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only

