

CM 23876 – POTASSIUM SORBATE 10% (10 ml/vl)

INTENDED USE

For the detection and isolation of acid resistant yeasts, *Zygosaccharomyces bailii* and *Zygosaccharomyces rouxii* in salads, sauces and dressings.

COMPOSITION

Ingredients	Concentration
Potassium sorbate	1g
Distilled water	10ml

(per vial sufficient for 200 ml medium)

INSTRUCTION FOR USE

Warm up refrigerated solution to room temperature. Aseptically add 1 ml to 1000 ml sterile, molten, cooled (45-50°C) APRY AGAR BASE along with 5 ml concentrated acetic acid. Mix well and pour into sterile Petri plate.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 Best Before	 QTY. Quantity	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 LOT/ B. NO. Lot / Batch Number	 Consults Instructions for Use	 QR Code	

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only