

CM 23878 – POTASSIUM TELLURITE 3.5% (1 ml/vl)

INTENDED USE

For selective isolation of Staphylococci and Corynebacteria.

COMPOSITION

Ingredients	Concentration
Potassium tellurite	0.350g
Distilled water	1ml

(per vial sufficient for 200ml medium)

INSTRUCTION FOR USE

Warm up the refrigerated contents of one vial to 45-50°C. Add aseptically 9 ml sterile distilled water, mix well and add 3 ml in 950 ml sterile, molten, cooled (45-50°C) BAIRD PARKER AGAR BASE / BAIRD PARKER AGAR BASE (VEG.) / BAIRD PARKER AGAR BASE (IS: 5887 (Part II) 1976, reaffirmed 2005.) / BAIRD PARKER AGAR BASE W/ SULPHA along with 50 ml Concentrated Egg Yolk Emulsion or 10 ml in 1000 ml GIOLITTI-CANTONI BROTH BASE / GIOLITTI-CANTONI BROTH BASE (ISO 6888-3:2003) / HOYLE MEDIUM BASE along with 50 ml of laked blood. Mix well and pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative Mallatier GmbH Buckhorn 10, 40143 Essen, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only

