

CM 23881 – RAPID LISTERIA SELECTIVE SUPPLEMENT

INTENDED USE

For selective enrichment of Listeria species from food samples.

COMPOSITION

Ingredients	Concentration
Amphotericin B	0.005g
Ceftazidime	0.010g
Nalidixic acid sodium salt	0.010g
Polymyxin B sulfate	38350 IU

(Per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial by adding of 4 ml of a 50:50 mixture of sterile demineralized water and ethanol (96 %). Mix gently until a homogeneous yellow suspension is obtained. Mix well and aseptically add to 475 ml of sterile, molten, cooled (45-50°C) LISTERIA ENRICHMENT BROTH BASE. Mix well and pour into sterile petri plates.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative Mueller GmbH Buckhorn 10, 48143 Münster, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only