

## CM 23883 – ROSOLIC ACID (0.1 g/vl)

### INTENDED USE

For isolation of coliform bacteria.

### COMPOSITION

| Ingredients  | Concentration |
|--------------|---------------|
| Rosolic acid | 0.100g        |

(per vial sufficient for 500 ml medium)

### INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 10 ml of 0.2 N sterile Sodium hydroxide. Mix well and aseptically add it to 990 ml of M-FC AGAR BASE / M-FC AGAR BASE (VEG.) / M-FC BROTH BASE before boiling. Add 10 ml of 1% Rosolic Acid to sterile, molten, cooled (45-50°C) M-FC AGAR BASE MODIFIED along with 50 mg Carbenicillin. Mix well and dispense as desired.

### STORAGE

Vials should be stored in sealed container at 2-8°C.

|                                                                                                                                         |                                                                                                                                                                                                       |                                                                                                            |                                                                                                                             |                                                                                                                           |                                                                                                       |
|-----------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------|
| <br>GMP<br>Good Manufacturing<br>Practices Certified | <br>IVD<br>For In Vitro Diagnostic Use                                                                             | <br>QTY.<br>Quantity    | <br>LOT/<br>B. NO.<br>Lot / Batch Number | <br>REF<br>Catalogue Number          | <br>Manufacturer |
| <br>Temperature Unit                                 | <br>EC REP<br>Authorized Representative<br><small>MedWet GmbH<br/>Buckenhof 16,<br/>48147 Münster, Germany</small> | <br>European Conformity | <br>QR<br>Code                           | <br>Consults Instructions<br>for Use | <br>Best Before  |

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

\*For Lab Use Only

