

CM 23887 – SALMONELLA SELECTIVE SUPPLEMENT

INTENDED USE

For isolation and identification of Salmonella spp.

COMPOSITION

Ingredients	Concentration
Novobiocin	2.5mg
Cefsulodin	6mg

(per vial sufficient for 500ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of 1 vial aseptically with 5ml of sterile distilled water. Mix gently to dissolve the contents completely. Aseptically add the rehydrated contents of one vial to 500ml of sterile, cooled (45-50°C) of CHROMOGENIC SALMONELLA AGAR. Mix well and dispense as desired.

STORAGE

Vial should be stored in sealed container at 2-8°C.

 Good Manufacturing Practices Certified	 For In Vitro Diagnostic Use	 Quantity	 Lot / Batch Number	 Catalogue Number	 Manufacturer
 Temperature Unit	 Authorized Representative MedWet GmbH Bockhorn 16, 48153 Bielefeld, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only