

CM 23898 – SULPHA SUPPLEMENT

INTENDED USE

For selective cultivation *Salmonella* species.

COMPOSITION

Ingredients	Concentration
Sodium sulphaacetamide	500mg
Sodium mandelate	125mg

(per vial sufficient for 500 ml medium)

INSTRUCTION FOR USE

Rehydrate the contents of one vial aseptically with 5.0 ml sterile distilled water. Mix well and aseptically add to 500 ml sterile, molten, cooled (45-50°C) BRILLIANT GREEN AGAR BASE, MODIFIED / BRILLIANT GREEN AGAR BASE W/ 1.2% AGAR / BRILLIANT GREEN AGAR BASE W/ PHOSPHATES / BRILLIANT GREEN AGAR BASE W/ PHOSPHATES (IS : 5887 (Part III) 1999, reaffirmed 2005) / TM 1512 - BRILLIANT GREEN AGAR W/ PHOSPHATES (ISO 6785:2001, ISO 6579:1981, ISO 3565:1975). Mix well and pour into sterile petri plates.

STORAGE

Vials should be stored in sealed container at 2-8°C.

 GMP Good Manufacturing Practices Certified	 IVD For In Vitro Diagnostic Use	 QTY. Quantity	 LOT/ B. NO. Lot / Batch Number	 REF Catalogue Number	 Manufacturer
 Temperature Unit	 EC REP Authorized Representative WuXi AppTec BioPharma Inc. 10113 Alameda, Germany	 European Conformity	 QR Code	 Consults Instructions for Use	 Best Before

NOTE: Please consult the Material Safety Data Sheet for information regarding hazards and safe handling Practices.

*For Lab Use Only